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
2012-06

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1998-06

GB - Instructions for use Your Peak Flow Meter

The éolys Peak Flow Meter is intended for use in the prevention and monitoring of Asthma by measuring the peak expiratory flow of adults and children.

Restrictions of use

Do not blow as you would do into a blowpipe; do not obstruct the movement of the slider; do not obstruct the air outlet; do not tighten or loosen the screws.

Instructions (see pictograms)

(1) Put the slider in the blue zone. (2) Breathe in as deeply as you can, with your chin up. (3) While biting gently on the mouthpiece. (4) Breathe out as HARD and as QUICKLY as you can. Repeat the operation twice, resting between each blow. (5) Remember to photocopy the grid before use. After having noted the date, record your best result of your 3 blows, indicated by the upper part of the cursor. Take a measurement in the morning (sun) and another one in the evening (moon). **IMPORTANT**: if you note a significant decrease in the results, contact your doctor.

Technical data

Accuracy +/- 10 L/min or +/-10 %. Resistance to the flow < 0,35 kPa//s. Readings are expressed as BTPS values.

Cleaning

Your appliance must be regularly cleaned. Immerse it in a slightly soapy tepid water (not warm and above all never boiling). Rinse and let it dry.

Maintenance

Your Peak Flow Meter does not require any particular maintenance apart from cleaning. Keep it in a clean and dry place.

Dysfunction and troubleshooting

You are advised to check your measuring appliance on a regular basis to make sure that it works properly. Check that the cursor can circulate normally and that all other parts remain fixed. If this is not the case you must send the appliance back to your supplier.

Disposal note

The device can be thrown away with the household waste.

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