



CareFusion

# SpiroUSB

Service manual




## Spiro USB - System Overview (Fig. 1)

The Spiro USB is a PC connected spirometer dedicated to work with Spirometry PC Software (SPCS).

It consists of a removable digital volume transducer (1) and a housing (2) containing a microprocessor control circuit and USB driver.

When testing a subject the transducer is inserted into the housing, which is plugged into a USB socket of a PC. The digital volume transducer is used to measure the subjects expired flow and volume in accordance with the operating manual.




## Transducer (Fig. 2)

The CareFusion digital volume transducer consists of an acrylic tube with a vane positioned between two swirl plates. The low inertia vane is attached to a stainless steel pivot which is free to rotate on two jewelled bearings mounted at the centre of the swirl plates. As air is passed through the transducer a vortex is created by the swirl plates which causes the vane to rotate in a direction dependant upon the direction of air flow. The number of rotations is proportional to the volume of air passed through the transducer and the frequency of rotation is proportional to the flow rate. The transducer housing consists of a main body which contains a pair of light emitting diodes (LED's) and phototransistors. The transducer is fixed to the mouthpiece holder which pushes into the main body and is captured by an "O" ring seal. The LED's produce infra red beams which are interrupted by the vane twice per revolution. This interruption is sensed by the phototransistors. The output from the collector of each phototransistor will be a square wave with a phase difference between the two of + or - 90 degrees depending upon the direction of flow.


There is no routine maintenance required for the transducer other than cleaning according to the instructions in the operating manual.

**CareFusion Digital Volume Transducer**



## Disassembly (Fig. 3)

1. Carefully remove both Spiro USB labels from the transducer housing.
2. Remove the screw under each label.
3. Pull apart the two halves of the housing as shown below:




## Reassembly (Fig. 4)

Please note: Do not use excessive force when reassembling.

1. Ensure the PCB is aligned correctly as shown in Fig. 4.
2. Pull the cable gland through the hole in housing A (already on the cable) and locate the centre moulding into the housing.
3. Line up the channel on the centre moulding with the screw hole of housing A.
4. When refitting housing B, locate the moulded bracket (with the nut) between the end of the channel and housing A and carefully pivot housing B around until the two housings meet (ensuring the blue LED fits into it's hole).
5. Replace the screws.
6. Carefully reposition the Spiro USB labels.

The PCB, LED'S and phototransistors will be accessible as shown below:

Fig. 4.



## **Microprocessor control circuit, See drawing 075-01 and 075-02**

The microprocessor control circuit monitors the transducer pulses, carries out the spirometry routines, and communicates with the PC via a USB driver under the control of it's internal program.

Power for the processor circuit is derived from the 5 volt USB power line. The power line is filtered by C12, L1 and C13. The filtered 5 volt is regulated down to 3.3 volts by the linear regulator, U7. U7 also provides a reset signal for the microprocessor.

The microprocessor, U3, is a Hitachi HD64F2318 16 bit microprocessor with 256K of flash memory and 8K of Ram. The system clock is supplied by 12MHz crystal, X1. There is also 512 Kbytes of external RAM, U8, used for storing pulses during a spirometry manoeuvre. The internal flash memory is used to store the microprocessor firmware.

Calibration data and system data is stored in an EEPROM, U2. Communication to the EEPROM is carried out using a two wire serial connection to pins 54 and 55 of the microprocessor. If the device is ever replaced, the unit will have to undergo factory recalibration.

Ambient temperature is monitored by a solid-state temperature sensor, U9. It communicates with the microprocessor via a one wire serial interface on pin 90.

The ambient temperature reading is used for adjusting inspiratory volume at ambient temperature to volume at body temperature.

The supply to the two series LED's, mounted inside the transducer housing, is provided through TR1, which is switched on by pin 4 of the microprocessor during a spirometry manoeuvre. Inside the transducer housing the two phototransistors used to detect the

interrupted infra-red beam are in open collector configuration. The pull up resistor for the two phototransistors is provided by R7 and R8.

Pulses from the phototransistor, TR2, are applied to the pulse timing input of the processor, pin 5, after being squared up by the action of the Schmitt inverter, U4. Pulses from the second phototransistor, TR3, after conditionings, U5, are applied to pin 6 of the microprocessor and are used to determine the direction of flow. The pulse count is used to determine the volume passed through the transducer since the start of the test and the pulse period is used to determine the flow at each volume increment.

The microprocessor communicates with the host PC via a USB interface, U6. U6 is connected to the microprocessor data bus and one address line, A0.

The 512k x 8 static RAM, U8, is located on the underside of the PCB and communicates with the processor on the 19-bit address bus and 8-bit data bus.


The speaker, J1, is directly connected to ports on pins 91 and 92 that are toggled at 1 KHz to generate the sound.

Parts List For: Spiro USB			
Drawing No. 075-00			Date: 21/8/2006
Revision No. 1.11			Page: 1 OF 2
Designation	Dwg./Manufacturers Part no.	JDE Part no.	Description
<b>36-ASS1041 PCB ASSEMBLY</b>			
U1	BU4S11	36-C1042	Single Schmitt NAND gate, SOT23-5 package
U2	24LC00-OT	36-C1025	128 bit serial EEPROM SOT23-5 package
U3	HD64F2318VTE25	36-C1057	Hitachi microcontroller
U4	INC-BU4SS584	36-C1043	Single Schmitt inverter, SOT23-5 package
U5	BU4S584	36-C1043	Single Schmitt inverter, SOT23-5 package
U6	INC-SL811HST-AXC	36-C1123	Cypress USB controller TQFP48 package
U7	MAX6349TL	36-C1103	3V3 regulator with integrated RESET control SOT23-5 package
U8	K6X4008T1F-VF70	36-C1065	Samsung 512k X 8 bit CMOS static RAM, SOL32/525 package
U9	DS18S20	36-C1049	Dallas semiconductor digital thermometer, T-O-92B package
U10	BU4S11	36-C1042	Single Schmitt NAND gate, SOT23-5 package
TR1	TRA-DTD113EK	36-TRA1045	Rohm NPN digital transistor SOT 23 package
TR2	BSS138	36-TRA1044	FET transistor, SOT-23 package
R1	100K0603	36-RES1044	100K surface mount 0.06 watt 1% resistor 0603 package
R2	100K0603	36-RES1044	100K surface mount 0.06 watt 1% resistor 0603 package
R3	1K0603	36-RES1086	1K surface mount 0.06 watt 1% resistor 0603 package
R4	10K0603	36-RES1051	10K surface mount 0.06 watt 1% resistor 0603 package
R5	10K0603	36-RES1051	10K surface mount 0.06 watt 1% resistor 0603 package
R6	1K0603	36-RES1086	1K surface mount 0.06 watt 1% resistor 0603 package
R7	15K0603	36-RES1212	15K surface mount 0.06 watt 1% resistor 0603 package
R8	15K0603	36-RES1212	15K surface mount 0.06 watt 1% resistor 0603 package
R9	22R0603	36-RES1210	22 Ohm surface mount 0.06 watt 1% resistor 0603 package
R10	22R0603	36-RES1210	22 Ohm surface mount 0.06 watt 1% resistor 0603 package
R11	150R0603	36-RES1071	1.5K surface mount 0.06 watt 1% resistor 0603 package
R12	120R0805	36-RES1064	120 Ohm surface mount 0.1 watt 1% resistor 0805 package
R13	100K0603	36-RES1044	100K surface mount 0.06 watt 1% resistor 0603 package
R14	100K0603	36-RES1044	100K surface mount 0.06 watt 1% resistor 0603 package
R15	100K0603	36-RES1044	100K surface mount 0.06 watt 1% resistor 0603 package
R16	100K0603	36-RES1044	100K surface mount 0.06 watt 1% resistor 0603 package
C1	33P0603	36-CAP1040	33pF 0603 package
C2	33P0603	36-CAP1040	33pF 0603 package
C3	100N0603	36-CAP1009	100nF 0603 package
C4	100N0603	36-CAP1009	100nF 0603 package
C5	100N0603	36-CAP1009	100nF 0603 package
C6	100N0603	36-CAP1009	100nF 0603 package
C7	1N0603	36-CAP1029	1nF ceramic 0603 package
C8	1N0603	36-CAP1029	1nF ceramic 0603 package
C9	100N0603	36-CAP1009	100nF ceramic 0603 package
C10	100N0603	36-CAP1009	100nF ceramic 0603 package
C11	4700N1206	36-CAP1047	4.7uF Ceramic 1206 package
C12	100N0603	36-CAP1009	100nF ceramic 0603 package
C13	1M0805	36-CAP1026	1uF AVX or Philips capacitor size 0805 X7R dielectric
L1	LQW18ANR22G00D	36-IND1023	Murata 220nH inductor, 0603 package
X1	CS1012.000MABJTR	36-XTL1002	Citizen 12MHz surface mount crystal, Digi-Key 300-8089-1-ND
CON1	IL-Z-5P-S125T3-E	36-CON1033	JAE 5 way pin header, Farnell 388-7170
CON2	53261-0490	36-CON1026	Molex 4 way pin header, 90 degree mount, Farnell 889-430
J1	PT-1550W	36-MEC1275	Dau Components piezo ceramic sounder
J2	05FKZ-SM1-1-TB	36-CON1067	JST ZIF FFC connector 5 way
LED2	SFH487-2-lead free	36-DIO1036	Osram infrared LED 3MM/T1
LED3	L934MBC	36-DIO1042	Kingbright 3mm blue LED, RS 247-1628
	075-03	36-PCB1014	PCB 075-03 ISS 3
<b>FLEXI CIRCUIT ASSEMBLY</b>			
	075-20	36-PCB1015	PCB 075-20 ISS 1.1
D1	SFH487-2-lead free	36-DIO1036	Osram infrared LED 3MM/T1
TR1	SFH309FA-6-lead free	36-TRA1006	Osram Photo transistor
TR2	SFH309FA-6-lead free	36-TRA1006	Osram Photo transistor



Parts List For: Spiro USB			
Drawing No. 075-00			Date: 21/8/2006
Revision No. 1.11			Page: 2 OF 2
Designation	Manufacturers Part No. (F=Farnell)	JDE Part no.	Description
	<b>039-18</b>	<b>36-TDX1043</b>	<b>TURBINE ASSEMBLY</b>
	031-17	36-MLD1388	2 off moulded turbine inner swirl plate
	031-18	36-MLD1389	2 off moulded turbine outer swirl plate
	039-12	36-MLD1431	Transducer Body - Reduced Cone Diameter
	039-02	36-MLD1344	Mouthpiece holder Grey
	039-14	36-MEC1168	Jewelled fixed bearing M3 thread
	039-16	36-MEC1168	Jewelled sprung bearing M3 thread
	039-17	36-ASS1152	27mm offset vane assembly
	039-15	36-MEC1228	27mm Pivot
	039-08	36-MEC1160	Flow deflector plate (039-08 iss 1.3)
	<b>075-14</b>	<b>36-CAB1016</b>	<b>CABLE ASSEMBLY</b>
	039-00	36-MLD1088	Transducer holder - body
	GUC-AO-1.8M	36-CAB1046	Glory Mark USB A Male to open Cable Assembly Grey
	075-05	36-MEC1186	Moulded cable gromett
	51021-0400	36-CON1020	Molex 4 way receptacle housing, Farnell 889-489
	50079-8100	36-CON1005	4 off Molex crimp terminal, 26/28 AWG female.
	CB01654	36-MEC1141	CPC 4.8mm black heatshrink sleeving as required
	150-152	36-FIX1073	Farnell tie wrap.
			<b>CASE COMPONENTS</b>
	039-00	36-MLD1088	Transducer holder - body
	075-04	36-MLD1294	Transducer body machining
	039-01	36-MLD1093	2 off transducer holder - LED housing
	BS028N170	36-MEC1394	BSL O ring 1.78 x ID 34.65mm Nitrile 70 Shore
	SCW-M2X4CSK	36-FIX1043	2 off M2 X 4mm countersunk posidrive screw
	NUT-M2	36-FIX1007	2 off M2 nut
	075-07	36-LAB1103	2 off Spiro USB logo label
	075-08	36-LAB1104	USB cable label
	PAD-2198232	36-ADH1017	3 off Sellotape double-sided self adhesive pad 12x25x1.5mm
			<b>SUNDRY ITEMS</b>
	GUC-AAF-3M	36-CAB1047	Glory Mark USB2.0 A Male to Female Cable Assembly Grey
	44B8008	36-LAB1179	CILS Wheelie bin self adhesive label 12.5 x 9.5mm





Title			Spiro USB - Circuit diagram, memory		
Size	Number			Revision	
A3	075-02			1.0	
Date:	9-Nov-2009		Sheet	of	
File:	C:\Documents and Settings\arthur.sadler\Desktop\Bsp075\075-02\Bsp075-02.dwg				

CareFusion UK 232 Ltd.,  
Quayside  
Chatham Maritime  
Kent ME4 4QY  
U.K.



JDE no. 36-MAN1152  
Drg no. 075-12  
Issue 1.1  
August 2009  
© CareFusion 2009

