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Le réflexe stapedial est mesuré at the static ear canal pressure qui produces the maximum membrane compliance, so reflex measurements are taken after the tympanogram is measured when the peak compliance pressure has been established. The reflex stimulus sound pressure level (SPL) is increased in steps until the change in compliance exceeds a predetermined threshold. The change in compliance at that level when the stimulus is applied is displayed as a plot against time.

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Lors du changement d'embout, assurz-vous que le joint est correctement inseré (posé) with its flat aligned with the flat on the probe tip. Push the probe tip over the boss and replace the nose cone. Make sure that the

Otowave 102 Operating Manual v.4

nose cone is screwed home firmly but do not over-tighten. Do not use any tools to tighten the nose cone.

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| MAIN MENU     | Ē      |
|---------------|--------|
| NEW TEST      |        |
| VIEW THE LAST | TEST   |
| DAILY CHECK   | Select |

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A typical tympanogram measurement and reflex test is carried out as follows.

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Present the ear tip to the ear and obtain a seal. If a good seal has been detected the following sequence of messages will be seen

| TESTING LEFT EAR      |
|-----------------------|
| "Equalising Pressure" |
|                       |
| Cancel                |
|                       |
| TESTING LEFT EAR      |
| "Pressure Settling"   |
|                       |
|                       |
| Cancel                |
|                       |
| TESTING LEFT EAR      |
| Seal Obtained 🗸       |
| Taking Tympanogram    |
| Cancel                |

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Page 14

Avant de commencer le test réflexe the ear canal pressure will be set to the value that gave the peak compliance during the tympanogram test. The instrument will then step through the tone frequencies and levels set in the CONFIGURATION menu searching for a reflex response:

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The display shows:

- The peak compliance, in ml (Pk)
- The pressure which gave the peak compliance in daPa
- The Gradient, in daPa (Gr)
- Le volume du canal auditif (ECV) en ml mesuré à 200 daPa.
- A plot of compliance against pressure.

Review the tympanogram to make sure that the peak compliance point selected by the Otowave is correct. If you are not satisfied you may select another peak using the  $\uparrow$  and  $\downarrow$  keys. The figures displayed will change to reflect the peak you select.

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If the reflex test was carried out the results will now be displayed:

Otowave 102 Operating Manual v.4

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- La fréquence de la mesure.
- "PASS" if a reflex was found, else "NR" (No Response).
- The level of the tone for which a reflex was first found.
- A plot of compliance against time.

If the reflex test was performed at more than one frequency use the  $\uparrow$  and  $\downarrow$  keys to view the results for the other frequencies.

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Otowave 102 Operating Manual v.4

After confirming that you wish to send data the message "Trying to connect..." will be displayed.

This will change to "Connection OK" when a connection to the computer has been made, and then to "Sending record…". A message will appear on the computer screen asking if you wish to accept the data. Click the "Yes to all" button and the data will be transferred.

When the data has been sent you will be returned to the previous menu.

Press  $\leftarrow$  at any time to stop sending data.

The transmitted results are placed in a folder called "Amplivox". By default this is placed on the current users desktop. If the folder already exists subsequent data will be saved in folders called "Copy 1 of Amplivox", "Copy 2 of Amplivox" and so on.

Each test is stored in a separate file within the folder. Files are named thus:

#### nnn\_DDMMYYYY\_HHMM.APX

where nnn is the identifier entered when you stored the test in the tympanometer (see section 6), or "xxx" if no identifier is available. DDMMYYYY is the date the measurement was saved and HHMM is the time the measurement was saved.

If you get a "Device not found" message while trying to send data check the following:

- The Otowave is pointing directly at the computers infra-red receiver from no more than 20cm away.
- The computer has its IrDA software properly installed and the interface enabled.
- If the computer has been in "Hibernate" mode the IrDA interface is not always re-enabled. Try restarting your computer.
- The IrDA adapter on your computer is compatible with the Otowave.

Turn the Otowave off and on again before trying to send the data again.

If communication is lost while sending the data you will get a "Link was unreliable" message. Press  $\leftarrow$  to cancel sending the data and start the operation again.

Page 22

If you see any other messages while sending data, turn the Otowave off and then on again. Try sending the data again. If the problem persists contact your Amplivox service centre.

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| data management                                   |  |
|---|--|
| LIST RECORDS                                      |  |
| DELETE RECORDS                                    |  |
| PRINT RECORDS<br>Back $\uparrow\downarrow$ Select |  |

Scroll down to see the remaining choice:

| data management                                     |
|---|
| DELETE RECORDS                                      |
| PRINT RECORDS                                       |
| SEND RECORDS TO PC Back $\uparrow\downarrow$ Select |

If you wish to work with the record of an individual test, select LIST RECORDS. All other options operate on groups of records.

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The sealing washer should be replaced if it shows signs of wear, or if a pressure leak is suspected.

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If the instrument is to be used at elevations above that specified in Section 2.1 re-calibration must be undertaken at the intended operating elevation.

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| Sub-menu              | Option     | Choix / Descriptions  |
|-----------------------|------------|---|
| NEW TEST              | SELECT EAR | Choisissez quelle(s) oreille(s) à<br>tester to test et commencer le test.<br>Un tympanogramme is taken<br>followed by reflex measurements, if<br>selected. On-screen messages &<br>LEDs indicate progress. Graphical<br>displays are shown automatically at<br>the end. |
| VIEW THE LAST<br>TEST | SELECT EAR | Recall the last stored test for the<br>selected ear. Shows the<br>tympanogram and reflex responses,<br>if available. Also allows the last test  |

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## 13. Messages d'erreur

| Message                                   | Meaning / Action                 |
|---|----------------------------------|
| PROBE NOT CLEAR                           | Examine the probe tip for        |
| Please ensure the probe is not blocked    | blockages. If necessary take it  |
| or obstructed                             | off and clean or replace it, see |
| PUMP ERROR.                               | section 5.3. If the problem      |
| Unknown pump fault. Restart the unit. If  | persists, contact your Amplivox  |
| problem persists, contact Amplivox        | service centre.                  |
| WARNING! CALIBRATION EXPIRED.             | The current date is later than   |
| Recalibration needed before further tests | the next calibration date. Check |
| are performed                             | that the clock is set to the     |

Page 30

|   | correct date. If so, arrange for<br>the instrument to be<br>recalibrated. Tests can still be<br>performed.   |
|---|--|
| "WARNING! BATTERIES LOW.<br>Remplacez les piles avant d'effectuer<br>des tests  | Remplacez les piles<br>immediatement, voir section<br>5.1  |
| Powering down   | L'Otowave est éteint because<br>the batteries are spent.<br>Replace the batteries.   |
| PUMP ERROR. Cannot determine pump<br>direction. If problem persists, contact<br>Amplivox                                    | Pump fault. If the fault persists contact your Amplivox service centre.  |
| PUMP ERROR. If problem persists, contact Amplivox   |  |
| Measurement timed out   | This occurs when the ear seal<br>check is set to THOROUGH if:<br>(i) The pump failed to achieve<br>the starting pressure within 4<br>seconds. This may be because<br>the probe was moved in the<br>ear.<br>(ii) The pressure failed to reach<br>-400 daPa within 12 seconds.<br>Retry the test. If the problem<br>persists, contact your Amplivox<br>service centre. |
| "WARNING! DEVICE UNCALIBRATED.<br>One or more default values require<br>recalibration before further tests are<br>performed | This message should never<br>normally be seen. If it persists<br>contact your Amplivox service<br>centre.  |
| WARNING! DEFAULTS RELOADED.<br>Default configuration settings reloaded.<br>Check before making new tests                    | Ce message ne devrait jamais<br>apparaître. Check all the<br>CONFIGURATION settings<br>before taking any<br>measurements. If the error<br>persists, contact your Amplivox<br>service centre.   |

Otowave 102 Operating Manual v.4

| ERROR                 | L'Otowave n'a pas pu envoyer     |
|-----------------------|----------------------------------|
| Transfer failed       | de donnée à l'ordinateur. Voir   |
| No device found       | section 8 pour plus de détails.  |
| ERROR                 |                                  |
| Transfer failed       |                                  |
| Link was unreliable   |                                  |
| WITHDRAW PROBE        | La sonde a bougé pendant la      |
|                       | prise mesure. Réinsérez la       |
|                       | sonde pour répéter le test       |
| Volume outside range  | The ear canal volume is above    |
| WITHDRAW PROBE        | the 5ml. This message also       |
|                       | occurs when the probe is not     |
|                       | properly inserted into the ear.  |
| Blocked probe         | The ear canal volume is below    |
| WITHDRAW PROBE        | 0.1ml. This message also         |
|                       | occurs when the probe tip is     |
|                       | blocked. Check that the probe    |
|                       | is correctly inserted into the   |
|                       | ear. Check that the probe is not |
|                       | blocked.                         |
| PLEASE REINSERT PROBE | The seal was lost. Reinsert the  |
|                       | probe to repeat the test.        |

## 14. Specifications techniques

| Tympanometrie                  |   |
|--------------------------------|---|
| Type d'instrument              | Meatus compensated tympanometer         |
| Analysis performed             | Compliance peak level (in ml); Pressure |
|                                | of same; Gradient (in daPa);            |
|                                | Ear Canal Volume (ECV) @ 200 daPa       |
| Probe tone levels and accuracy | 226Hz +/- 2%; 85dB SPL +/-2dB over      |
|                                | range 0.2ml to 5ml                      |
| Pressure levels and accuracy   | +200daPa to -400daPa +/-10daPa or       |
|                                | +/-10% (whichever is larger) over range |
| Ear volume measurement range   | 0.1ml to 5ml +/- 0.1ml or +/-5%         |
| and accuracy                   | (whichever is larger) over entire range |
| Sweep speed                    | Typically 200-300daPa/sec; dependant    |

## 14.1. Performance

Page 32

|  | on ear/cavity volume   |  |
|--|--|--|
| Pressure limits (safety cutout)              | +600 to -800 daPa  |  |
| Number of samples stored                     | 100 per tympanogram  |  |
| Reflex measurements                          |  |  |
| Modes de prise de mesure                     | Ipsilateral  |  |
| Reflex tone levels and accuracy              | 102-1: 1000Hz +/-2%, configurable over<br>range 60dB to100dB HL +/-2dB,<br>referenced to 2ml calibration volume;<br>Compensates for measured ear volume  |  |
|  | 102-4: 500Hz, 1kHz, 2kHz, 4kHz (4 kHz restricted to 95dBHL)  |  |
| Reflex measurement range and accuracy        | 0.01ml to 0.5ml +/-0.01ml configurable<br>in 0.01ml steps  |  |
| Number of reflex levels                      | 4, in 5 dB or 10 dB steps  |  |
| Analyse du réflexe                           | Reflex pass/fail at each level tested;<br>maximum amplitude of each reflex<br>(seen on printed report & PC report);<br>pressure at which reflex was performed  |  |
| Pressure used for reflex measurement         | Pressure at Tympanogram peak, or 0<br>daPa (Always and Prompt Before Each<br>Test modes)   |  |
| Reflex level cut-off                         | Optionally, Auto-stop when reflex found  |  |
| Reflex threshold detection                   | Configurable 0.01 – 0.50 ml in 0.01 ml increments  |  |
| Reflex tone duration                         | 0.6 seconds  |  |
| Number of records stored in Patient Database | 30   |  |
| Data storage                                 | Any recording can be stored once the tympanogram is viewed. Patient Initials (A-Z, 0-9, "-") must be entered before storage.   |  |
| Data held                                    | Patient Initials, Tympanogram and<br>Reflex graphs and analysis for Left Ear<br>and/or Right Ear, Time and Date of<br>recording, which ears were tested,<br>whether or not the record has been<br>printed and/or sent to a PC, parameters<br>used for analysis, 128 bit Globally<br>Unique Identifier (GUID) |  |

Otowave 102 Operating Manual v.4

|  | 1  |
|--|--|
| Display mode                               | Records listed in reverse chronological<br>order (latest first), with indication of data<br>stored as described above  |
| Real Time Clock                            |  |
| Time stamps                                | Time and date stamp applied to all recordings, and to the last calibration date  |
| Backup power supply                        | > 30 days without main batteries fitted  |
| Printing                                   |  |
| Supported printer                          | MCP8830  |
| Interface                                  | Infra-red, IrDA hardware, 9600 baud  |
| Information imprimée                       | Space for patient & clinician's details,<br>Tympanogram analysis parameters,<br>Tympanogram, Reflex analysis<br>parameters, Reflex graph, Serial<br>Number of device, Last and Next Due<br>Calibration dates |
| Serial Interface to PC                     |  |
| Interface                                  | OBEX (Object Exchange) service<br>running on top of IrDA stack. Auto-<br>selects rate between 9600 - 115200<br>baud.   |
| Serial Interface to PC                     |  |
| Information envoyée                        | Patient header, full left and right ear data.  |
| Power Supply                               |  |
| Battery Types                              | <ul> <li>4 Alkaline AA Cells or;</li> <li>4 NiMH rechargeable</li> <li>NiMH batteries which must be of greater than 2.3 Ah capacity.</li> </ul>  |
| Warm-up period                             | None at room temperature   |
| Number of recordings from one set of cells | Approx 200 (Alkaline) or 100 (NiMH)  |
| Auto power-off delay                       | 90 or 180 S  |
| Idle current                               | 70mA   |
| Current while testing                      | 230mA  |
| Physical                                   |  |
| Display                                    | 128 x 64 pixels / 8 lines of 21 characters   |

Page 34

| Dimensions                     | 190mm long x 80mm wide x 40mm high |
|--------------------------------|------------------------------------|
|                                | excluding probe                    |
|                                | 210mm long including probe         |
| Poids (sans les piles)         | 285 g                              |
| Poids (avec piles)             | 380 g                              |
| Environmental                  |                                    |
| Operating temperature range    | +10 C to +40 C                     |
| Operating humidity range       | 10% – 90% RH, non-condensing       |
| Operating atmospheric pressure | 980 – 1040 mb                      |
| range                          |                                    |
| Storage temperature range      | - 40 C to +70 C                    |
| Storage humidity range         | 10% – 90% RH, non-condensing       |
| Storage atmospheric pressure   | 900 – 1100 mb                      |
| range                          |                                    |
| Standards conformance          |                                    |
| Safety                         | EN60601-1: 1990                    |
| EMC                            | EN60601-1-2: 2001                  |
| Performance                    | EN60645-5: 2005, Type 2            |
|                                | Tympanometer                       |
| CE mark                        | To the Medical Device Directive    |

## 14.2. Equipment classification

Type de protection contre les chock Degree of protection against electric shock Degree of protection against ingress of water Mode of operation Equipment mobility Internally Powered Type BF applied part Not protected Continuous operation Portable

The Otowave 102 Tympanometer is classified as a Class IIa device under Annex IX (Section 1) of the Medical Devices Directive. It is intended for transient use as a screening tympanometer instrument.

## 14.3. Symboles

Definition: Type BF equipment – equipment providing a particular degree of protection against shock, particularly regarding allowable LEAKAGE current and reliability of the protective earth connection (if present).

Otowave 102 Operating Manual v.4

Définition: Attention, consult accompanying documents.

Page 36

## **Commander consommables et Accessoires**

To order consumables, additional accessories and to replace detachable parts that have been damaged, please contact Amplivox for current prices and delivery charges. The items available are listed below:

| Stock No. | Description                                   |
|-----------|---|
| T517      | Probe tip                                     |
| T518      | Sealing Washer                                |
| T018      | Test Chamber Volume, 0.5 and 1 ml             |
| T10       | Ear Tip Set                                   |
| T101      | Embout oriculaire Otowave 6mm                 |
| T102      | Embout oriculaire Otowave 7mm                 |
| T103      | Embout oriculaire Otowave 8mm                 |
| T104      | Embout oriculaire Otowave 9.5mm               |
| T105      | Embout oriculaire Otowave 11mm                |
| T106      | Embout oriculaire Otowave 12.5mm              |
| T107      | Embout oriculaire Otowave 14mm                |
| T108      | Embout oriculaire Otowave 16mm                |
| T109      | Embout oriculaire Otowave 18mm                |
| B132      | Malette de transport                          |
| MANOW     | Manuel d'utilisation de l'Amplivox Otowave    |
| A091      | Imprimante thermique                          |
| C01       | Thermal Printer rolls (20metres)              |
| T91       | ACTiSYS Infrared USB adapter                  |
| T003      | Amplivox NOAH Impedance module                |
| T004      | Amplivox NOAH Impedance module + IrDA PC port |

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Page 38

# Guidance and manufacturer's declaration – electromagnetic immunity (1)

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| Immunity test                      | IEC 60601<br>test level            | Compliance<br>level | Electromagnetic<br>environment –<br>guidance  |
|------------------------------------|------------------------------------|---------------------|---|
| Electrostatic<br>Discharge (ESD)   | ±6 kV contact                      | ±6 kV contact       | Floors should be<br>wood, concrete or   |
| IEC 61000-4-2                      | TO KV All                          |                     | are covered with<br>synthetic material,<br>the relative humidity<br>should be at least<br>30% |
| Electrical fast<br>transient/burst | ±2 kV for<br>power supply<br>lines | Not applicable      | Not applicable  |
| IEC 61000-4-4                      |                                    |                     |   |
|                                    | ±1 kV for<br>input/output<br>lines |                     |   |
| Surge                              | ±1 kV<br>differential              | Not applicable      | Not applicable  |
| IEC 61000-4-5                      | mode                               |                     |   |
|                                    | ±2 kV                              |                     |   |
|                                    | common                             |                     |   |
|                                    | moue                               |                     |   |

Otowave 102 Operating Manual v.4

| Voltage dips,<br>short<br>interruptions and<br>voltage variations<br>on power supply<br>input lines<br>IEC 61000-4-11 | <5% U <sub>T</sub><br>(>95% dip in<br>U <sub>T</sub> ) for 0.5<br>cycle<br>40% U <sub>T</sub><br>(60% dip in<br>U <sub>T</sub> ) for 5<br>cycles<br>70% U <sub>T</sub><br>(30% dip in<br>U <sub>T</sub> ) for 25<br>cycles<br><5% U <sub>T</sub><br>(>95% dip in<br>U <sub>T</sub> ) for 5 sec | Not applicable | Not applicable   |
|---|--|----------------|--|
| Power frequency<br>(50/60 Hz)<br>magnetic field<br>IEC 61000-4-8  | 3 A/m  | 3 A/m          | Power frequency<br>magnetic fields<br>should be at levels<br>characteristic of a<br>typical location in a<br>typical commercial<br>or hospital<br>environment. |
| $N \cup I = \cup_T$ is the a.c. mains voltage prior to the application of the test level                              |  |                |  |

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| Immunity<br>test        | IEC 60601<br>test level | Compliance<br>level       | Milieu électromagnétique – guide d'information   |
|-------------------------|-------------------------|---------------------------|--|
|                         |                         |                           | Portable and mobile RF<br>communications equipment<br>should be used no closer to<br>any part of the Otowave 102<br>Tympanometer, including<br>cables, than the<br>recommended separation<br>distance calculated from the<br>equation applicable to the<br>frequency of the transmitter. |
| Radiated<br>RF          | 3 V/m<br>80MHz to       | 3 V/m                     | Recommended separation distance  |
| IEC 2.5GHz<br>61000-4-3 |                         | d = 1.2√P 80MHz to 800MHz |  |
|                         |                         |                           | d = 2.3√P 800MHz to 2.5GHz   |
|                         |                         |                           | where P is the maximum<br>output power rating of the<br>transmitter in Watts (W)<br>according to the transmitter<br>manufacturer and d is the<br>recommended separation<br>distance in metres (m).   |
|                         |                         |                           | Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each   |

Otowave 102 Operating Manual v.4

| Guide d'information et déclaration du fabricant – immunité<br>électromagnétique (2) |   |  |   |
|---|---|--|---|
|   |   |  | frequency range. <sup>b</sup>   |
|   |   |  | Interference may occur in the vicinity of equipment marked with the following symbol: |
|   |   |  | (((●)))   |
| NOTE '  | 1 At 80MHz and 800N   | MHz, the higher f                      | requency range applies.   |
| NOTE 2<br>propaga<br>and peo  | 2 These guidelines m<br>ation is affected by at<br>ople.  | ay not apply in a<br>osorption and ref | Il situations. Electromagnetic<br>lection from structures, objects                    |
| а   | Les champs de force from fixed transmitters, such as base stations<br>for radio (cellular/cordless) telephones and land mobile radios,<br>amateur radio, AM and FM radio broadcast and TV broadcast<br>cannot be predicted theoretically with accuracy. To assess the<br>electromagnetic environment due to fixed RF transmitters, an<br>electromagnetic site survey should be considered. If the measured<br>field strength in the location in which the Otowave 102<br>Tympanometer is used exceeds the applicable RF compliance<br>level above, the Otowave 102 Tympanometer should be observed<br>to verify normal operation. If abnormal performance is observed,<br>additional measures may be necessary, such as re-orienting or<br>relocating the Otowave 102 Tympanometer. |  |   |
| b   | over the frequency range in the should be less than a   | ange 150 kHz to<br>3 V/m.              | 80 MHz, field strengths   |

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| maximum output p                                   |   | incations equipment | it.            |  |
|--|---|---------------------|----------------|--|
| Rated<br>maximum<br>output power<br>of transmitter | Separation distance according to frequency of transmitter |                     |                |  |
|  | 150 kHz to 80   | 80 MHz to 800       | 800 MHz to 2.5 |  |
| W  | MHz   | MHz                 | GHz            |  |
|  | d = 1.2√P   | d = 1.2√P           | d = 2.3√P      |  |
| 0.01   | 0.12  | 0.12                | 0.23           |  |
| 0.1  | 0.38  | 0.38                | 0.73           |  |
| 1  | 1.2   | 1.2                 | 2.3            |  |
| 10   | 3.8   | 3.8                 | 7.3            |  |
| 100  | 12  | 12                  | 23             |  |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

**<u>NOTE 1</u>** At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Otowave 102 Operating Manual v.4