

# PistonXP Release notes

Piston Kft.

---

The logo consists of the letters 'RN' in a bold, white, sans-serif font, centered within a solid black square.

**Version: 2.1-rc.11**

**Release date: 2025-06-05**

# TABLE OF CONTENTS

<b>Updating PistonXP</b>	<b>3</b>
Updating from versions before 2.10.0.0698 . . . . .	3
Updating from versions before 1.84.0.0519 . . . . .	3
Updating from versions before 1.82.0.0458 . . . . .	3
<b>Changelog</b>	<b>4</b>
[2.10.0.0709] - 2025-06-05 . . . . .	4
Fixed . . . . .	4
Changed . . . . .	5
Added . . . . .	5
[2.00.0.0680] - 2024-01-02 . . . . .	6
Fixed . . . . .	6
Changed . . . . .	6
Added . . . . .	7
Removed . . . . .	7
[1.84.0.0594] - 2020-10-01 . . . . .	7
Fixed . . . . .	8
[1.84.0.0579] - 2020-05-04 . . . . .	8
Fixed . . . . .	8
Added . . . . .	8

# UPDATING PISTONXP



**Before installing the update, make a backup copy of the PistonXP installation directory!**



**For OEM versions:**

**Do not use the **Online Update** button in **PistonXP** **Setup** **Maintenance**! This would download and update to the latest official version, with Piston logo.**

There are several ways to update PistonXP:

- Run Setup ▶ xpx\_update.exe from the installer DVD and click the Install button.

, or

- Autorun the installer DVD (runs Autorun.exe)
- Select **Software** ▶ **Install PistonXP and drivers** (runs Setup ▶ xpx\_fullsetup.exe)
- On the **Select Components** screen, choose **Update the software only** and click the **Install** button

, or you can simply just reinstall it to the same directory.

After the installer has finished the update, run PistonXP.exe and depending on the version you are updating from, do the following extra steps:

## Updating from versions before 2.10.0.0698

- **Audiometer:** Connect your audiometer, then set **PistonXP** ▶ **Setup** ▶ **Audiometer** ▶ **Standard** ▶ **Read from device** and **PistonXP** ▶ **Setup** ▶ **Audiometer** ▶ **Headset** ▶ **Read from device**.

## Updating from versions before 1.84.0.0519

- **CPET:** Perform a sensor characteristics measurement: **PistonXP** ▶ **Setup** ▶ **CPET** ▶ **Measurement of gas sensor characteristics** ▶ **Start calibration**.

## Updating from versions before 1.82.0.0458

- Create admin password.
- Create user and setup auto login if needed.

# CHANGELOG

All notable changes to this project will be documented in this file. The format is based on [Keep a changelog](#)

## [2.10.0.0709] - 2025-06-05

---

### Fixed

- CPET window: FVC data source indication - Ticket 205
- Prediction method indication for stored parameters - Ticket 242
- Double confirmation of “store data” at exit - Ticket 226
- Base spirometry: An extra Pred column is visible - Ticket 211
- Remove empty “results” tab for others test types than FVC - Ticket 212
- Exact prediction algorithm type not stored when Custom is selected - Ticket 220
- Trend Analysis bug fixes - Ticket 222
- Missing parameters on CPET report - Ticket 192
- Parameter display options are language dependent in INI file - Ticket 217
- List only used prediction algorithms on report - Ticket 219
- Mutex issue - Ticket 310
- Table headers on printed reports of CPET diagnoses - Ticket 312
- Displaying import/export menu item - Ticket 314
- Stability fixes on ILLA interface for Eolys - Ticket 315
- Stability and accuracy fixes on SPIRO interface for Eolys - Tickets 317, 322
- Stability fixes for plethysmograph tests - Ticket 321
- Stability fixes for CPET - Ticket 324
- Stability fixes for SpO2 and trend analysis - Ticket 325
- Character set compatibility issue - Ticket 326
- Minor refinements and fixes on provocation test - Ticket 328
- Unresponsive OK button of Settings window on PCs without COM ports - Ticket 331
- Usability improvements of Comment editor - Ticket 338
- Fixes on importing far eastern characters via GDT - Ticket 340
- Bug fix for Bluetooth heart rate monitors - Ticket 346
- Stability fixes for PinkFlow zeroing - Ticket 351
- Stability fixes for CPET - Ticket 348
- Date filter fix on patient form’s detailed search - Ticket 359
- CSV format fix on SPIRO interface for Eolys - Ticket 360, 364
- Stability fixes, refinements and customization options for flow calibration - Ticket 363
- Unification of device serial number on test reports - Ticket 370
- Minor stability fixes - Ticket 372
- FEV1/IVC and FEV1/VC calculation fix - Ticket 379
- Motorized pump plug-in issue fix - Ticket 375

- Predicted value mark-up, calculation and re-calculation fixes - Ticket 380
- Minor bug fixes on FOT communication services
- Fix for copy-paste feature on parameter visibility matrix - Ticket 385
- Fix for hangs while generating reports with unset default printer - Ticket 388
- Fix for MVV parameter calculation - Ticket 398
- Fix for patient filtering in PostgreSQL databases - Ticket 399
- Fix for maximal respiratory pressure gauge - Ticket 400
- PostgreSQL backup fix - Ticket 412

### Changed

- Low screen resolution warning - Ticket 255
- Grading display moved to parameters table - Ticket 228
- Only notify to calibrate if any of the devices is set to do so - Ticket 208
- Data import/export for MediSAP, Eolys - Ticket 318
- GUI and translation refinements for CPET - Ticket 324
- Customized OEM installers and DB preset for Eolys - Ticket 333
- Minor code optimizations - Ticket 373
- GDT export is changed to multi-body (multi-sentence) to improve compatibility - Ticket 343
- Refining calculation of end-tidal gas concentrations - Ticket 396
- GDT parameter ordering - Ticket 407

### Added

- Set parameter visibility at once - Ticket 256
- Added warnings to Measurement Details Form - Ticket 233
- Configuration file formatted automatically - Ticket 223
- Guidelines on Wasserman panels - Ticket 259
- Audio Threshold table "Police Nationale", eolys - Ticket 246
- Explanatory table for FVC grading system - Ticket 229
- Predefined composite prediction algorithms - Ticket 227
- NHANES III Predicted values - Ticket 174
- PADOA HIS interface for Eolys - Ticket 199
- Basic spirometry: auto stop when no flow detected - Ticket 210
- More details are added for calibration step descriptions - Ticket 209
- Additional transcutaneous pCO<sub>2</sub>/pO<sub>2</sub> monitor support - Ticket 313
- Limited PostgreSQL 16 support - Ticket 317
- Error messages during HIS (GDT) issues are detailed - Ticket 327
- Option for customization of CPET warm-up parameters by protocol - Ticket 334
- New features in Comment editor, new measurement count badges on Main menu - Ticket 352
- Vyaire/Padoa command line interface for Eolys - Ticket 347, 409
- Usability and ergonomics enhancements - Ticket 356
- INI file backup - Ticket 309
- AMEDTEC ECGpro integration - Ticket 329
- Support for truncated patient ID while loading back earlier CPET via GDT - Ticket 358

- New spirometry features (recording SpO2 and HR, animation, zero count-down, Z-Score on report) - Ticket 361
- GDT compatibility support for Innomed CardioPC/E 3.2 - Ticket 362
- Support for misformatted CSV input on SPIRO interface for Eolys - Ticket 367
- FOT mode switch warning added - Ticket 371
- INI file validator - Ticket 330
- Adapted connectivity for Bluetooth LE devices - Ticket 355
- Customizable semi-automatic PRE/POST function selection for Spirometry - Ticket 376
- Distinguish inactive doctors (closed praxis), change patient-doctor association - Ticket 377
- Customizable PEmax/Plmax on report, new and enhanced HL7 export features for SOTE - Ticket 381
- Out-of-canvas drawing highlighted for Spirometry - Ticket 382
- Patient merge feature - Ticket 383
- GUI improvements for Jinbo Medical - Ticket 389
- Enriched CPET data export for data evaluation - Ticket 393
- FOT device support update and intrabreath mode - Ticket 394, 405
- Pulse interface on spirometry for Eolys - Ticket 395
- Manual upload of predicted spirometry parameters on CPET for SSEM - Ticket 397
- Protocol customization via GUI on CPET for TSO Medical - Ticket 403
- Customizable noise filtering on PinkFlow - Ticket 411
- Added support for GSS CO2 and Hummingbird O2 sensors - Ticket 408

## **[2.00.0.0680] - 2024-01-02**

---

### **Fixed**

- Calibration time lag calculation bug - Ticket 60
- DataID, MasterID anomaly in CHDATA table - Ticket 93
- PWG software couldn't connect - Ticket 125
- Plethy - Diff measurement screen switch bug - Ticket 15
- DLCO test sometime switch to IntraBreath mode - Ticket 76
- PistonXP freezes with patients height = 60cm - Ticket 94
- Threadmill protocol steps table bug - Ticket 176
- Measurement order mixed up - Ticket 183
- Comment window disappears - Ticket 126
- Gaw, sGaw parameter units - Ticket 150
- PXP update database path from ini - Ticket 27
- CPET report 670 message - Ticket 118
- Audiometer USB removal - Ticket 190
- Minor bugfixes

### **Changed**

- Side menu redesign - Ticket 184
- Minimal screen resolution is now 1600x900
- Mouth pressure measurement enhancement - Ticket 12
- Lower PinkFlow low flow threshold (+1% accuracy) - Ticket 181

- Restrictions on CPET TCP interface - Ticket 142
- Graph zoom with CTRL + mouse wheel
- Graph autozoom with double click
- FOT CV instead of coherence - Ticket 244
- Lactate diagram standardization - Ticket 275

### Added

- ATS/ERS 2019 Compliance - Ticket 77
- GLI 2020 and 2021 predictions update - Ticket 186
- New **ETSPLs for Piston PDD-401 circumaural audiometric headphones**
- Audiometer ETSPL and PELTOR correction table read from device - Ticket 180
- Serial number and software version on report - Ticket 8
- DB update tool - Ticket 70
- Writing comments for patients - Ticket 121
- Tlco Hb, Dlco Hb and Kco Hb parameters - Ticket 187
- Tlco Hb, Dlco Hb and Kco Hb **predictions**
- FRC parameter in TGV measurement - Ticket 137
- Diffusion capacity test VA parameter - Ticket 150
- LLN and Z-Score diagram on report - Ticket 14
- SpO2 test: show previous data - Ticket 107
- Printing comments for calibration report - Ticket 127
- Russian language support - Ticket 128
- VE/VCO2 slope on the result tab on Slopes section - Ticket 116
- Increase decrease load from keyboard - Ticket 117
- Automatic flow calibration (PCS-1500/M) - Ticket 182
- Adaptive, scriptable CPET protocol framework - Ticket 197
- Continuous save during CPET examination - Ticket 225
- FOT adjustable measurement duration - Ticket 244
- FOT calibration tube impedance calculation - Ticket 244
- FOT fres, Ax4, Ax5 predictions - Ticket 299
- CPET dead space calculation (VD/VT-Jones and VD/VT ABG) - Ticket 274
- Digital signature

### Removed

- Legacy FTDI device compatibility removed
- Legacy MCHPUSB device compatibility removed
- DBSetup removed (compatibility issues)
- User manuals removed from the installer (still available on the installation CD)

## [1.84.0.0594] - 2020-10-01

---

Released **only for CPET** system! Known bugs with Plethysmograph and FOT.  
This release is only compatible with **BTL CardioPoint version 2.33**.

---

**Fixed**

- Rhinomanometer graph and PRE/POST – Ticket 114
- PO2 and PCO2 parameter calculation – Ticket 115
- HRR doesn't correspond with HR predicted - HR max – Ticket 108
- OUES calculation - wrong value – Ticket 110

**[1.84.0.0579] - 2020-05-04**

---

Released **only for CPET** system! Known bugs with Plethysmograph.

**Fixed**

- Comments disappearing totally – Ticket 51
- Error opening record by GDT – Ticket 61
- CardioPoint incomplete measurement report error – Ticket 4
- Warm-up error message pops up strangely – Ticket 52
- Weird record – Ticket 53
- Detecting device disconnection – Ticket 59
- Comments printing – Ticket 3
- CPET settings store - Ticket 16
- Average values displayed incorrectly – Ticket 40
- Piston report does not merge to CardioPoint report automatically - Ticket 41
- RCP threshold doesn't place to intersection of two lines - - Ticket 42
- CPET: find in DB – Ticket 48
- English for Translation repetition – Ticket 54
- PXP merge doesn't work when updating 498 version – Ticket 56
- Incorrect icon of Piston appears – Ticket 57
- CPET report contains BMR value, makes no sense – Ticket 58
- Calibration time lag calculation – Ticket 60
- Store port setting – Ticket 62
- CPET - Tab Table is not refreshed after the test is done – Ticket 66
- Table items saved in template change the position – Ticket 71
- Results table (finished test) in Table card must be refreshed to see the values – Ticket 72
- VO2max, VCO2max options not stored – Ticket 73
- New patient not selected (PSQL) – Ticket 74
- Starting second instance will crash – Ticket 46

**Added**

- Added some contents to the User's Manual – according to TÜV requests – Ticket 49