CareFusion

SpiroUSB

Operating Manual – English

Operating Manual – English - USA

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Handbok – Svenska

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Contents

I. Introduction - English	13
Package contents	14
Contraindications	15
Warnings and Cautions	15
Indication for Use	16
Installing SPCS	17
Installing SpiroUSB device driver	17
Compatibility status	17
Driver installation	17
Operation	18
Looking after your SpiroUSB Spirometer	18
Calibration Check	19
Cleaning Instructions	20
External Surfaces of the Spirometer	20
Cleaning accessories	20
Cleaning the Transducer	22
Servicing	23
Product Lifetime	23
Trouble Shooting Information	24
Safety Designation per IEC 60601-1	25
Electromagnetic Compatibility (EMC) to EN60601-1-2:2007	27
Symbols	32
Specification of the SpiroUSB	33
Spirometry Measurements:	34
Consumables/ Accessories	36
Customer contact information	38
II. Introduction – English - USA	42
Package contents	43

Contraindications	43
Warnings and Cautions	44
Indication for Use.....	45
Installing SPCS	46
Installing SpiroUSB device driver	46
Compatibility status.....	46
Driver installation.....	46
Operation.....	47
Looking after your SpiroUSB Spirometer	47
Calibration Check	48
Cleaning Instructions	49
External Surfaces of the Spirometer	49
Cleaning accessories	49
Cleaning the Transducer.....	51
Servicing	52
Product Lifetime.....	52
Trouble Shooting Information.....	53
Safety Designation per IEC 60601-1	54
Electromagnetic Compatibility (EMC) to EN60601-1-2:2007.....	56
Symbols	61
Specification of the SpiroUSB	62
Spirometry Measurements:.....	63
Consumables/ Accessories	65
Customer contact information	67
III. Introduction - Français	70
Contenu de l'emballage	71
Contre-indications	72
Avertissements et mises en garde	72
Indications d'utilisation	74

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Operating manual

I. Introduction - English

The SpiroUSB is a PC connected spirometer dedicated to work with the comprehensive Spirometry PC Software (SPCS).

SPCS is a fully Windows™ compatible spirometry system that interfaces seamlessly with the SpiroUSB providing many display options and includes powerful reporting and database facilities.

SPCS has many advanced features including display of real time respiratory traces, predicted values, patient database, and the ability to carry out pre and post bronchodilator and post steroid testing.


The SpiroUSB uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. This transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

Package contents

The SpiroUSB is packaged in a sturdy carrying case containing this manual and the following items (Fig.1):

1. Carefusion Digital Volume Transducer
2. SpiroUSB transducer housing.

Together with SPCS, extension cable, disposable cardboard mouthpieces and nose clip.

Contraindications

- Acute disorders affecting test performance (e.g Vomiting , nausea, vertigo)
- Recent eye surgery (increases in intraocular pressure during spirometry)
- Oral or facial pain exacerbated by a mouthpiece
- recent myocardial infarction
- post- operative thoracic surgery patients
- hyperventilation syndrome

Note: Extensive exhalation might lead to syncope

Warnings and Cautions

The following terms are used as follows in this manual

CAUTION: Possibility of injury or serious damage

WARNING: conditions or practices that could result in personal injury

Note: Important information for avoiding damage to the instrument or facilitating operation of the instrument.

Note: Patients below the age of 4 may struggle to perform spirometry correctly and reproducibly

Note: The device should only be used by trained and qualified personnel.

CAUTION: Read the manual before use

WARNING: The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anaesthetic mixtures or in oxygen rich environments.

CAUTION: Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement."

CAUTION: Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement."

PLEASE NOTE: The product you have purchased should not be disposed of as unsorted waste. Please utilise your local WEEE collection facilities for the disposal of this product.

PLEASE NOTE: Degree of protection against Ingress of Water is IPX0

WARNING: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth

CAUTION: When you connect the SpiroUSB to other equipment, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the SpiroUSB only to printers and computers that comply with IEC 60601-1 / ANSI/AAMI ES60601-1:2005 / CAN/CSA-C22.2 No. 60601-1:14.

Indication for Use

The SpiroUSB spirometer is intended, for prescription use only, to measure the maximal volume and flow of air that can be moved in and out of a patient's lungs. The system is intended for use with paediatric (4 to 17 years) and adult (18 to 99 years) patients in hospitals, physician offices, laboratories and occupational health testing environments.

Installing SPCS

SPCS is supplied on CD ROM with your SpiroUSB. Install SPCS according to the instructions in the SPCS quick start guide.

Installing SpiroUSB device driver

The driver enables CareFusion Ltd USB devices to work with USB enabled PC software on Windows based machines.

Compatibility status

There are two versions of the USB driver. The driver used depends on the Operating System you are using. There is a driver for Windows 32 bit Operating Systems and for 64 bit.

The following operating systems are supported:

- Windows 7 32 & 64 bit Operating Systems
- Windows 8 32 & 64 bit Operating Systems

The driver will not be supported under older Windows Operating Systems.

Driver installation

If software is installed that provides support for USB, the drivers may be pre-installed, allowing automatic installation of the drivers when the device is connected to the USB port. If the device does not install automatically, download the driver and install the device manually. The driver can be downloaded from the following homepage:

www.carefusion.co.uk/our-products/respiratory-care/cardio-pulmonary-diagnostics/pulmonary-function-testing/spirometers/spirometry-software-and-firmware-downloads

Operation

Simply connect the SpiroUSB to a spare USB port on your PC with the supplied SPCS installed. SPCS will automatically detect that the SpiroUSB has been connected to a USB port and will run, ready for you to start your spirometry session. There is no need to remove the SpiroUSB between sessions.

Please note that with certain PC configurations, even if a compatible operating system is being used, SPCS may not be launched automatically when the SpiroUSB is plugged in. In this case SPCS may be run in the usual way using the Windows on-screen Start button.

The SpiroUSB transducer is supplied with a two-meter cable. If this length is insufficient then a remote USB hub should be used.

Remote USB hubs are readily available from computer accessory suppliers. For some recent computers the USB signal will work over a longer distance using an extension cable. An extension cable has been supplied for your convenience, if required. If communication problems are encountered using the extension cable, then a remote USB hub must be used.

Note: Keep the PC and Monitor out of reach of the patient at all times.

Looking after your SpiroUSB Spirometer

Please observe the following precautions:

- Avoid exposing the SpiroUSB to direct sunlight during use.
- Avoid operating the spirometer in dusty conditions or near to heating appliances or radiators.
- Do not keep the spirometer in a damp place or expose it to extremes of temperature.
- Do not direct the transducer holder towards a strong light source whilst operating the spirometer.

Calibration Check

The spirometer is calibrated to read in litres at body temperature, Barometric Pressure Saturated with water vapour (BTPS).

The calibration should remain stable indefinitely, unless the transducer is physically damaged, and the unit should not require re-calibration. However, to ensure the correct functioning of the unit we do recommend that a calibration check is performed after the transducer was removed for cleaning.

Pressing the 'Verify Calibration' button in SPCS invokes the start of a device search and the option Calibration Check. The device search effectively scans the PC's comms/usb ports for any attached device. On finding a device the screen will change to the main verify calibration test screen. The option to perform an expiratory only and expiratory and inspiratory calibration check, and a single volume calibration or more advanced calibration check at three separate flow rates, can be chosen in the calibration screen.

Connect a 3L syringe to the transducer with the minimum of adapters and empty it by pushing the handle fully in.

Note: It is recommended that the transducer is disinfected prior to a calibration check or a MicroGard® filter is used during the procedure.

[illegible]

Cleaning Instructions

Disinfection of contaminated parts is only effective after having them carefully pre-cleaned. CareFusion recommend the tested solution of PeraSafe sterilizing powder (36-SSC5000A) for pre-cleaning and disinfection. If a different solution is used, please follow the given manufacturer's instructions.

The device must not be exposed to solvents e.g. alcohol, chloride

CAUTION: Unplug your SpiroUSB before cleaning.

External Surfaces of the Spirometer

CAUTION: Do not attempt to wash or immerse the SpiroUSB transducer housing in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

The external housing of the spirometer may be wiped with sterile wipes or a damp cloth that has been immersed in a cold sterilising solution when required.

Cleaning accessories

With the use of a filter MicroGard® IIB (V-892381), MicroGard® IIC (V-892384) or a One-way Safety mouthpiece (36-SST1250) for each patient, cleaning for the components in patient's gas path is recommended once a month.

When using the disposable cardboard mouthpiece (adult: 36-PSA1000, pediatric: 36-PSA1200) without a filter and under the prerequisite that the patient was instructed only to exhale into the transducer, the following parts have to be cleaned once a day: adult/ pediatric adapter, transducer.

With any other use as described above all contaminated parts must be disinfected between patients.


Important note: Used single patient nose clip pads, mouthpieces, MicroGard® IIB filter and MicroGard® IIC filter must be disposed off immediately after the use.

If there are changes on the material surfaces (cracks, brittleness) the respective parts must be disposed off.

Cleaning the Transducer

The transducer requires no routine maintenance or servicing. However, if you wish to sterilise or clean the transducer it may be removed by means of the following procedure:

1. Rotating the turbine transducer anti-clockwise until the locating pip lines up with the small rectangular cut-out in the housing as shown below.
2. Gently pull the transducer away from the housing.
3. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in a cold sterilising solution e.g. Perasafe (Cat No. 36-SSC5000A) for a maximum of 10 minutes (Alcohol and chloride solutions should be avoided).
4. After cleaning/sterilising, the transducer should be rinsed briefly in distilled water and dried.
5. Re-assemble the mouthpiece holder.

CAUTION: Do not attempt to wash or immerse the SpiroUSB transducer housing in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

Servicing

A full service manual including circuit diagrams and parts list is available upon request from the manufacturer.

If your unit requires service or repair please see page 38 for contact details.

Product Lifetime

The SpiroUSB spirometer is designed for a product lifetime of 5 years.

Trouble Shooting Information

Should you encounter problems operating your Spiro USB meter, please consult the table below:

Problem	Possible Cause	Solution
Does not register a blow	Head assembly or cable broken	Replacement of head assembly or return the unit for servicing
Blows are inverted on the display	Head assembly or cable broken	Replacement of head assembly or return the unit for servicing
Blows tracing ends abruptly although patient still exhaling	Turbine sticking	Clean turbine in warm soapy water or sterilising solution, if problem continues a replacement turbine may be required
Calibration procedure failed or can not be completed	Turbine may be faulty	Repeat calibration procedure, if problem still persists, replace turbine or return unit for servicing
	Turbine not fitted tightly to calibration syringe	Ensure the syringe is fitted to the turbine using the adapter supplied part # 861427. (supplied with the syringe)
	Calibration syringe does not have an inspiratory seal or seal is leaking.	Ensure you are using manufacturer recommended syringe.
	Shaft of the syringe being pushed down	The syringe should be emptied and filled with one smooth stroke, avoid pushing down on the shaft or banging at the end of each manoeuver.

Safety Designation per IEC 60601-1

Type of protection against electrical shock	Powered by computers that comply with IEC 60601-1 / ANSI/AAMI ES60601-1:2005 / CAN/CSA-C22.2 No. 60601-1:14.
Degree of protection against electrical shock	Type B applied part
Power Equipment	Powered by computers that comply with IEC 60601-1 / ANSI/AAMI ES60601-1:2005 / CAN/CSA-C22.2 No. 60601-1:14.
Degree of Electrical connection between equipment and Patient	Equipment designed as non-electrical connection to the patient.
Degree of mobility	Transportable
Mode of operation	Continuous

Classifications according to IEC 60601-1

SpiroUSB	Applied part, type B
Volume Transducer	Applied part, type B

WARNING: No modification of this equipment is allowed.

CAUTION: When you connect other equipment to the unit, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the SpiroUSB only to printers and computers that comply with IEC 60601-1 / ANSI/AAMI ES60601-1:2005 / CAN/CSA-C22.2 No. 60601-1:14.

WARNING: The user must not touch any voltage-carrying parts and the patient at the same time or the operator must not create a “bridge” between the device I/O ports and the patient by simultaneously touching both.

WARNING: Do not connect devices that are not specified as part of the system.

WARNING: To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.

Electromagnetic Compatibility (EMC) to EN60601-1-2:2007

WARNING: use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation.

The SpiroUSB has been tested to EN60601-1-2:2007, regarding it's ability to operate in an environment containing other electrical/electronic equipment (including other medical devices).

The purpose of this testing is to ensure that the SpiroUSB is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of the SpiroUSB.

Despite the testing of the SpiroUSB that has been undertaken, normal operation of the SpiroUSB can be affected by other electrical/electronic equipment and portable and mobile RF communications equipment.

As the SpiroUSB is medical equipment, special precautions are needed regarding EMC (electromagnetic compatibility).

It is important that the SpiroUSB is configured and installed/put into service, in accordance with the instructions/guidance provided herein and is used only in the configuration as supplied.

Changes or modifications to the SpiroUSB may result in increased emissions or decreased immunity of the SpiroUSB in relation to EMC performance.

The SpiroUSB should be used only with the USB extension cable provided by Carefusion, which is referenced in the accessories section of this manual. This cable should not be used with devices other than the SpiroUSB. If the cable is extended by the user using non approved cables, this may result in an increased level of emissions or decreased level of immunity, in relation to the SpiroUSB's EMC. Use of the USB extension cable with devices other than the SpiroUSB, may result in an increased level of emissions or decreased level of immunity, in relation to the other devices's EMC.

The SpiroUSB has an essential performance – the product should continue to operate correctly. In the unlikely event of a Fast Transient / ESD event occurring, the device should be reset and located away from the source of interference.

WARNING: The SpiroUSB should not be used adjacent to or stacked with other equipment. If adjacent or stacked use with other equipment is necessary, the SpiroUSB and the other equipment should be observed/monitored, to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The SpiroUSB is intended for use in the electromagnetic environment specified below. The customer or the user of the SpiroUSB should assure that it is used in such an environment

Emission Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The SpiroUSB uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Group B	The SpiroUSB is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC61000-3-2	Not Applicable (USB powered)	
Voltage fluctuations / flicker emissions IEC61000-3-3	Not Applicable (USB powered)	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The SpiroUSB is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or the user of the SpiroUSB should ensure that it is used in such an environment

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment -Guidance
Electrostatic discharge (ESD) IEC61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	Supply line tests not applicable ± 1 kV for USB lead	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not Applicable (USB Powered)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short	<5% U_T (>95 % dip in U_T)	Not Applicable (USB)	Mains power quality should be that of a typical

interruptions and voltage variations on power supply input lines IEC61000-4-11	For 0.5 cycle 40% U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) For 5 s	Powered)	commercial or hospital environment. If the user of the SpiroUSB requires continued operation during power mains interruptions, it is recommended that the SpiroUSB be powered from an uninterruptable power supply or a battery.
Power frequency (50/60Hz) Magnetic field IEC61000-4-8	3 A/m	3 A/m	If incorrect operation occurs, it may be necessary to position the SpiroUSB further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The SpiroUSB is intended for use in the electromagnetic environment specified below. The customer or the user of the SpiroUSB should assure that it is used in such an environment

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the SpiroUSB, including any cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance (d) $d = 1.2\sqrt{P}$
Radiated RF	3 V/m	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz

IEC61000-4-3	80 MHz to 2.5 GHz		$d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SpiroUSB is used exceeds the applicable RF compliance level above, the SpiroUSB should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the SpiroUSB			
b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the SpiroUSB

The SpiroUSB is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the SpiroUSB can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SpiroUSB as recommended below, according to the maximum output power of the communications equipment

Rated Maximum Output Power of Transmitter in Watts (W)	Separation Distance in Meters (m) according to Frequency of Transmitter		
	150 KHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people

Symbols

Type B device

In accordance with Directive 93/42/EEC

0086

Disposal in compliance with WEEE

Consult the instructions for use

Manufacturer

Date of Manufacture

Serial number

Direct current

Single patient use

Rx only

Federal U.S. law restricts this device to sale by or on the order of a physician.

Specification of the SpiroUSB

General	
Power supply:	Input 5V 0.1A
Dimensions:	50x60x90mm
Weight:	130g
Operating Temperature:	0 to +40° C
Operating Humidity:	30% to 90% RH
Operating Pressure:	650 to 1060 hPa
Storage & Transport temperature:	-20 to + 70° C
Storage & Transport humidity:	10% to 90% RH
Storage & Transport Pressure:	650 to 1060 hPa
Predicted Values:	Various depends on National preference
Transducer:	Carefusion Bi-Directional Digital Volume.
Resolution:	10ml volume 0.03l/s flow
Accuracy:	+/-3%. To ATS recommendations – Standardisation of Spirometry 1994 update for flows and volumes.
Measurements:	See SPCS operating manual for details of all the spirometry measurements available.

Spirometry Measurements:

Relaxed Expiratory Vital Capacity (VC)

Forced Expired Volume in 0.75 seconds (FEV_{0.75})

Forced Expired Volume in 1 second (FEV₁)

Forced Expired Volume in 3 second (FEV₃)

Forced Expired Volume in 6 seconds (FEV₆)

Forced Vital Capacity (FVC)

Peak Expiratory Flow Rate (PEF)

FEV_{0.75} as a percentage of VC (FEV_{0.75}/VC)

FEV_{0.75} as a percentage of FVC (FEV_{0.75}/FVC)

FEV₁ as a percentage of VC (FEV₁/VC)

FEV₁ as a percentage of FVC (FEV₁/FVC)

FEV₃ as a percentage of VC (FEV₃/VC)

FEV₃ as a percentage of FVC (FEV₃/FVC)

FEV_{0.75} as a percentage of FEV₆ (FEV_{0.75}/FEV₆)

FEV₁ as a percentage of FEV₆ (FEV₁/FEV₆)

Maximum Expired Flow at 75% of FVC remaining (MEF75)

Maximum Expired Flow at 50% of FVC remaining (MEF50)

Maximum Expired Flow at 25% of FVC remaining (MEF25)

Mean Mid-Expiratory Flow Rate (MMEF)

Forced expiratory flow at 50% of volume as a percentage of VC (FEF50/VC)

Forced expiratory flow at 50% of volume as a percentage of FVC (FEF50/FVC)

Maximal voluntary ventilation indicated (MVV_(ind))

Forced inspired volume in 1 second (FIV₁)

Forced inspiratory Vital Capacity (FIVC)

Peak Inspiratory Flow Rate (PIF)

FIV₁ as a percentage of FIVC (FIV₁/FIVC)

Forced inspiratory flow at 25% of inhaled volume (FIF25)

Forced inspiratory flow at 50% of inhaled volume (FIF50)

Forced inspiratory flow at 75% of inhaled volume (FIF75)

Forced expiratory flow at 50% of volume as a percentage of FIF50 (FEF50/FIF50)

The time taken between 25% and 75% of the forced expired volume (MET2575)

Forced Expiratory Time (FET)

Tidal Volume (TV)

Expiratory reserve volume (ERV)

Inspiratory reserve volume (IRV)

Inspiratory capacity (IC)

Expiratory Relaxed vital capacity (EVC)

Inspiratory vital capacity (IVC)

Breathing frequency rate (FR)

Inspiratory time (Ti)

Expiratory time (Te)

Ti as a % of total breath time (Ti/Ttot)

Tidal volume as a % of Ti (TV/Ti)

Breath Rate BR

Breathing Time B.T

Volume Tidal VT

Expiratory Time - average time of expiration per breath in seconds Te

Inspiratory Time - average time of inspiration per breath in seconds Ti

Total Tidal Breath Time in Seconds TTOT = Ti + Te

Ratio of Average Expiratory and Inspiratory Breaths Ti/Te

Average Time of Expiration per Breath as a ratio to the Total Tidal Breath Time Ti/TTOT

Tests per subject: 5 VC manoeuvre

8 FVC manoeuvres

Predicted Values: Various – depends upon national preference

Transducer: CareFusion Bi-Directional Digital Volume.

Consumables/ Accessories

Cat. No.	Description
36-PSA1000	Adult Disposable Mouthpieces (500 per box)
36-SST1000	One-way Safety Mouthpieces (500 per box)
36-SST1250	One-way Safety Mouthpieces (250 per box)
36-PSA1200	Paediatric Disposable Mouthpieces (250 per box)
36-PSA1100	Paediatric Adaptor
V-892381	MicroGard® IIB with integrated mouthpiece (50 per box)
V-892384	MicroGard® IIC (50 per box)
36-SSC5000A	PeraSafe Sterilising Powder 81g (to make up 5 litres of solution)
36-VOL2104	Nose Clips (pack of 5)
892121	Nose clip pads "foam material", 100 pcs.
36-SM2125	3 Litre Calibration Syringe
V-861449	Silicone adapter "oval" for MicroGard® IIB
861427	Calibration syringe adapter
36-SPC1000	Spirometry PC Software
36-CAB1047	USB Extension cable
36-ASS1141	SpiroUSB desktop stand
36-TDX1043	Turbine Transducer

For further information or to place an order for Disposables/
Supporting products please contact CAREFUSION, your local
distributor or view our website www.carefusion.com/micromedical

PLEASE NOTE: USE ONLY CAREFUSION ACCESSORIES

Please Note: Information in this manual is subject to change without notice and does not represent commitment on the part of CareFusion 232 UK Ltd. The software may be used or copied only in accordance with the terms of that agreement. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of CareFusion 232 UK Ltd.

Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551

Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.internationalsales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: support.rt.eu@carefusion.com

SpiroUSB

Operating manual

II. Introduction – English - USA

The SpiroUSB is a PC connected spirometer dedicated to work with the comprehensive Spirometry PC Software (SPCS).

SPCS is a fully Windows™ compatible spirometry system that interfaces seamlessly with the SpiroUSB providing many display options and includes powerful reporting and database facilities.

SPCS has many advanced features including display of real time respiratory traces, predicted values, patient database, and the ability to carry out pre and post bronchodilator and post steroid testing.


The SpiroUSB uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. This transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

Package contents

The SpiroUSB is packaged in a sturdy carrying case containing this manual and the following items (Fig.1):

1. Carefusion Digital Volume Transducer
2. SpiroUSB transducer housing.

Together with SPCS, extension cable, disposable cardboard mouthpieces and nose clip.

Contraindications

- Acute disorders affecting test performance (e.g. Vomiting , nausea, vertigo)
- Recent eye surgery (increases in intraocular pressure during spirometry)
- Oral or facial pain exacerbated by a mouthpiece
- recent myocardial infarction
- post- operative thoracic surgery patients
- hyperventilation syndrome

Note: Extensive exhalation might lead to syncope

Warnings and Cautions

The following terms are used as follows in this manual

CAUTION: Possibility of injury or serious damage

WARNING: conditions or practices that could result in personal injury

Note: Important information for avoiding damage to the instrument or facilitating operation of the instrument.

Note: Patients below the age of 4 may struggle to perform spirometry correctly and reproducibly

Note: The device should only be used by trained and qualified personnel.

CAUTION: Read the manual before use

WARNING: The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anesthetic mixtures or in oxygen rich environments.

CAUTION: Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement."

CAUTION: Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement."

PLEASE NOTE: The product you have purchased should not be disposed of as unsorted waste. Please utilize your local WEEE collection facilities for the disposal of this product.

PLEASE NOTE: Degree of protection against Ingress of Water is IPX0

WARNING: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth

CAUTION: When you connect the SpiroUSB to other equipment, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the SpiroUSB only to printers and computers that comply with IEC 60601-1 / ANSI/AAMI ES60601-1:2005 / CAN/CSA-C22.2 No. 60601-1:14.

Indication for Use

The SpiroUSB spirometer is intended, for prescription use only, to measure the maximal volume and flow of air that can be moved in and out of a patient's lungs. The system is intended for use with pediatric (4 to 17 years) and adult (18 to 99 years) patients in hospitals, physician offices, laboratories and occupational health testing environments.

Installing SPCS

SPCS is supplied on CD ROM with your SpiroUSB. Install SPCS according to the instructions in the SPCS quick start guide.

Installing SpiroUSB device driver

The driver enables CareFusion Ltd USB devices to work with USB enabled PC software on Windows based machines.

Compatibility status

There are two versions of the USB driver. The driver used depends on the Operating System you are using. There is a driver for Windows 32 bit Operating Systems and for 64 bit.

The following operating systems are supported:

- Windows 7 32 & 64 bit Operating Systems
- Windows 8 32 & 64 bit Operating Systems

The driver will not be supported under older Windows Operating Systems.

Driver installation

If software is installed that provides support for USB, the drivers may be pre-installed, allowing automatic installation of the drivers when the device is connected to the USB port. If the device does not install automatically, download the driver and install the device manually. The driver can be downloaded from the following homepage:

www.carefusion.co.uk/our-products/respiratory-care/cardio-pulmonary-diagnostics/pulmonary-function-testing/spirometers/spirometry-software-and-firmware-downloads

Operation

Simply connect the SpiroUSB to a spare USB port on your PC with the supplied SPCS installed. SPCS will automatically detect that the SpiroUSB has been connected to a USB port and will run, ready for you to start your spirometry session. There is no need to remove the SpiroUSB between sessions.

Please note that with certain PC configurations, even if a compatible operating system is being used, SPCS may not be launched automatically when the SpiroUSB is plugged in. In this case SPCS may be run in the usual way using the Windows on-screen Start button.

The SpiroUSB transducer is supplied with a two-meter cable. If this length is insufficient then a remote USB hub should be used. Remote USB hubs are readily available from computer accessory suppliers. For some recent computers the USB signal will work over a longer distance using an extension cable. An extension cable has been supplied for your convenience, if required. If communication problems are encountered using the extension cable, then a remote USB hub must be used.

Note: Keep the PC and Monitor out of reach of the patient at all times.

Looking after your SpiroUSB Spirometer

Please observe the following precautions:

- Avoid exposing the SpiroUSB to direct sunlight during use.
- Avoid operating the spirometer in dusty conditions or near to heating appliances or radiators.
- Do not keep the spirometer in a damp place or expose it to extremes of temperature.
- Do not direct the transducer holder towards a strong light source whilst operating the spirometer.

Calibration Check

The spirometer is calibrated to read in liters at body temperature, Barometric Pressure Saturated with water vapour (BTPS).

The calibration should remain stable indefinitely, unless the transducer is physically damaged, and the unit should not require re-calibration. However, to ensure the correct functioning of the unit we do recommend that a calibration check is performed after the transducer was removed for cleaning.

Pressing the 'Verify Calibration' button in SPCS invokes the start of a device search and the option Calibration Check. The device search effectively scans the PC's commas/usb ports for any attached device. On finding a device the screen will change to the main verify calibration test screen. The option to perform an expiratory only and expiratory and inspiratory calibration check, and a single volume calibration or more advanced calibration check at three separate flow rates, can be chosen in the calibration screen.

Connect a 3L syringe to the transducer with the minimum of adapters and empty it by pushing the handle fully in.

Note: It is recommended that the transducer is disinfected prior to a calibration check or a MicroGard® filter is used during the procedure.

Patient
There is currently no patient selected

File Data Calibrate Help

Verify Calibration

Enter User Name:
Unit Name:
Device Serial Number:
Syringe Size: 3 Litres
Select Manoeuvre: ☐ Expiratory Only ☐ Expiratory and Inspiratory
Select Calibration Verification Requirements: ☐ Simple ☐ AT-SERS

Flow (L/s) vs Volume (L) graph showing a red line at 0.

Stage	Status	IDEAL	MEAS
Band (row 0 - 9 9H4)			
Test 1	EXP	Waiting	

View History Start Test Exit Verify Calibration Next Band

Cleaning Instructions

Disinfection of contaminated parts is only effective after having them carefully pre-cleaned. CareFusion recommend the tested solution of PeraSafe sterilizing powder (36-SSC5000A) for pre-cleaning and disinfection. If a different solution is used, please follow the given manufacturer's instructions.

The device must not be exposed to solvents e.g. alcohol, chloride

CAUTION: Unplug your SpiroUSB before cleaning.

External Surfaces of the Spirometer

CAUTION: Do not attempt to wash or immerse the SpiroUSB transducer housing in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

The external housing of the spirometer may be wiped with sterile wipes or a damp cloth that has been immersed in a cold sterilizing solution when required.

Cleaning accessories

With the use of a filter MicroGard® IIB (V-892381), MicroGard® IIC (V-892384) or a One-way Safety mouthpiece (36-SST1250) for each patient, cleaning for the components in patient's gas path is recommended once a month.

When using the disposable cardboard mouthpiece (adult: 36-PSA1000, pediatric: 36-PSA1200) without a filter and under the prerequisite that the patient was instructed only to exhale into the transducer, the following parts have to be cleaned once a day: adult/ pediatric adapter, transducer.

With any other use as described above all contaminated parts must be disinfected between patients.


Important note: Used single patient nose clip pads, mouthpieces, MicroGard® IIB filter and MicroGard® IIC filter must be disposed off immediately after the use.

If there are changes on the material surfaces (cracks, brittleness) the respective parts must be disposed off.

Cleaning the Transducer

The transducer requires no routine maintenance or servicing. However, if you wish to sterilize or clean the transducer it may be removed by means of the following procedure:

1. Rotating the turbine transducer anti-clockwise until the locating pip lines up with the small rectangular cut-out in the housing as shown below.
2. Gently pull the transducer away from the housing.
3. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in a cold sterilizing solution e.g. Perasafe (Cat No. 36-SSC5000A) for a maximum of 10 minutes (Alcohol and chloride solutions should be avoided).
4. After cleaning/sterilizing, the transducer should be rinsed briefly in distilled water and dried.
5. Re-assemble the mouthpiece holder.

CAUTION: Do not attempt to wash or immerse the SpiroUSB transducer housing in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

Servicing

A full service manual including circuit diagrams and parts list is available upon request from the manufacturer.

If your unit requires service or repair please see page 67 for contact details.

Product Lifetime

The SpiroUSB spirometer is designed for a product lifetime of 5 years.

Trouble Shooting Information

Should you encounter problems operating your Spiro USB meter, please consult the table below:

Problem	Possible Cause	Solution
Does not register a blow	Head assembly or cable broken	Replacement of head assembly or return the unit for servicing
Blows are inverted on the display	Head assembly or cable broken	Replacement of head assembly or return the unit for servicing
Blows tracing ends abruptly although patient still exhaling	Turbine sticking	Clean turbine in warm soapy water or sterilizing solution, if problem continues a replacement turbine may be required
Calibration procedure failed or cannot be completed	Turbine may be faulty	Repeat calibration procedure, if problem still persists, replace turbine or return unit for servicing
	Turbine not fitted tightly to calibration syringe	Ensure the syringe is fitted to the turbine using the adapter supplied part # 861427. (supplied with the syringe)
	Calibration syringe does not have an inspiratory seal or seal is leaking.	Ensure you are using manufacturer recommended syringe.
	Shaft of the syringe being pushed down	The syringe should be emptied and filled with one smooth stroke, avoid pushing down on the shaft or banging at the end of each maneuver.

Safety Designation per IEC 60601-1

Type of protection against electrical shock

Powered by computers that comply with IEC 60601-1 / ANSI/AAMI ES60601-1:2005 / CAN/CSA-C22.2 No. 60601-1:14.

Degree of protection against electrical shock

Type B applied part

Power Equipment

Powered by computers that comply with IEC 60601-1 / ANSI/AAMI ES60601-1:2005 / CAN/CSA-C22.2 No. 60601-1:14.

Degree of Electrical connection between equipment and Patient

Equipment designed as non-electrical connection to the patient.

Degree of mobility

Transportable

Mode of operation

Continuous

Classifications according to IEC 60601-1

SpiroUSB

Applied part, type B

Volume Transducer

Applied part, type B

WARNING: No modification of this equipment is allowed.

CAUTION: When you connect other equipment to the unit, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the SpiroUSB only to printers and computers that comply with IEC 60601-1 / ANSI/AAMI ES60601-1:2005 / CAN/CSA-C22.2 No. 60601-1:14.

WARNING: The user must not touch any voltage-carrying parts and the patient at the same time or the operator must not create a “bridge” between the device I/O ports and the patient by simultaneously touching both.

WARNING: Do not connect devices that are not specified as part of the system.

WARNING: To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.

Electromagnetic Compatibility (EMC) to EN60601-1-2:2007

WARNING: use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation.

The SpiroUSB has been tested to EN60601-1-2:2007, regarding its ability to operate in an environment containing other electrical/electronic equipment (including other medical devices).

The purpose of this testing is to ensure that the SpiroUSB is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of the SpiroUSB.

Despite the testing of the SpiroUSB that has been undertaken, normal operation of the SpiroUSB can be affected by other electrical/electronic equipment and portable and mobile RF communications equipment.

As the SpiroUSB is medical equipment, special precautions are needed regarding EMC (electromagnetic compatibility).

It is important that the SpiroUSB is configured and installed/put into service, in accordance with the instructions/guidance provided herein and is used only in the configuration as supplied.

Changes or modifications to the SpiroUSB may result in increased emissions or decreased immunity of the SpiroUSB in relation to EMC performance.

The SpiroUSB should be used only with the USB extension cable provided by Carefusion, which is referenced in the accessories section of this manual. This cable should not be used with devices other than the SpiroUSB. If the cable is extended by the user using non approved cables, this may result in an increased level of emissions or decreased level of immunity, in relation to the SpiroUSB's EMC. Use of the USB extension cable with devices other than the SpiroUSB, may result in an increased level of emissions or decreased level of immunity, in relation to the other devices' EMC.

The SpiroUSB has an essential performance – the product should continue to operate correctly. In the unlikely event of a Fast Transient / ESD event occurring, the device should be reset and located away from the source of interference.

WARNING: The SpiroUSB should not be used adjacent to or stacked with other equipment. If adjacent or stacked use with other equipment is necessary, the SpiroUSB and the other equipment should be observed/ monitored, to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The SpiroUSB is intended for use in the electromagnetic environment specified below. The customer or the user of the SpiroUSB should assure that it is used in such an environment

Emission Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The SpiroUSB uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Group B	The SpiroUSB is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC61000-3-2	Not Applicable (USB powered)	
Voltage fluctuations / flicker emissions IEC61000-3-3	Not Applicable (USB powered)	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The SpiroUSB is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or the user of the SpiroUSB should ensure that it is used in such an environment

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment -Guidance
Electrostatic discharge (ESD) IEC61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	Supply line tests not applicable ± 1 kV for USB lead	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not Applicable (USB Powered)	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5% U_T (>95 % dip in U_T) For 0.5 cycle 40% U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) For 5 s	Not Applicable (USB Powered)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the SpiroUSB requires continued operation during power mains interruptions, it is recommended that the SpiroUSB be powered from an uninterruptable power supply or a battery.
Power frequency (50/60Hz) Magnetic field IEC61000-4-8	3 A/m	3 A/m	If incorrect operation occurs, it may be necessary to position the SpiroUSB further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The SpiroUSB is intended for use in the electromagnetic environment specified below. The customer or the user of the SpiroUSB should assure that it is used in such an environment

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the SpiroUSB, including any cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance (d) $d = 1.2\sqrt{P}$

Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a . Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SpiroUSB is used exceeds the applicable RF compliance level above, the SpiroUSB should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the SpiroUSB			
^b . Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the SpiroUSB

The SpiroUSB is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the SpiroUSB can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SpiroUSB as recommended below, according to the maximum output power of the communications equipment

Rated Maximum Output Power of Transmitter in Watts (W)	Separation Distance in Meters (m) according to Frequency of Transmitter		
	150 KHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people

Symbols

Type B device

In accordance with Directive 93/42/EEC

0086

Disposal in compliance with WEEE

Consult the instructions for use

Manufacturer

Date of Manufacture

Serial number

Direct current

Single patient use

Rx only

Federal U.S. law restricts this device to sale by or on the order of a physician.

Specification of the SpiroUSB

General	
Power supply:	Input 5V 0.1A
Dimensions:	50x60x90mm
Weight:	130g
Operating Temperature:	0 to +40° C
Operating Humidity:	30% to 90% RH
Operating Pressure:	650 – 1060 hPa
Storage & Transport temperature:	-20 to + 70° C
Storage & Transport humidity:	10% to 90% RH
Storage & Transport Pressure:	650 – 1060 hPa
Predicted Values:	Various depends on National preference
Transducer:	Carefusion Bi-Directional Digital Volume.
Resolution:	10ml volume 0.03l/s flow
Accuracy:	+/-3%. To ATS recommendations – Standardization of Spirometry 1994 update for flows and volumes.
Measurements:	See SPCS operating manual for details of all the spirometry measurements available.

Spirometry Measurements:

Relaxed Expiratory Vital Capacity (VC)

Forced Expired Volume in 0.75 seconds (FEV_{0.75})

Forced Expired Volume in 1 second (FEV₁)

Forced Expired Volume in 3 second (FEV₃)

Forced Expired Volume in 6 seconds (FEV₆)

Forced Vital Capacity (FVC)

Peak Expiratory Flow Rate (PEF)

FEV_{0.75} as a percentage of VC (FEV_{0.75}/VC)

FEV_{0.75} as a percentage of FVC (FEV_{0.75}/FVC)

FEV₁ as a percentage of VC (FEV₁/VC)

FEV₁ as a percentage of FVC (FEV₁/FVC)

FEV₃ as a percentage of VC (FEV₃/VC)

FEV₃ as a percentage of FVC (FEV₃/FVC)

FEV_{0.75} as a percentage of FEV₆ (FEV_{0.75}/FEV₆)

FEV₁ as a percentage of FEV₆ (FEV₁/FEV₆)

Maximum Expired Flow at 75% of FVC remaining (MEF₇₅)

Maximum Expired Flow at 50% of FVC remaining (MEF₅₀)

Maximum Expired Flow at 25% of FVC remaining (MEF₂₅)

Mean Mid-Expiratory Flow Rate (MMEF)

Forced expiratory flow at 50% of volume as a percentage of VC (FEF₅₀/VC)

Forced expiratory flow at 50% of volume as a percentage of FVC (FEF₅₀/FVC)

Maximal voluntary ventilation indicated (MVV_(Ind))

Forced inspired volume in 1 second (FIV₁)

Forced inspiratory Vital Capacity (FIVC)

Peak Inspiratory Flow Rate (PIF)

FIV₁ as a percentage of FIVC (FIV₁/FIVC)

Forced inspiratory flow at 25% of inhaled volume (FIF₂₅)

Forced inspiratory flow at 50% of inhaled volume (FIF₅₀)

Forced inspiratory flow at 75% of inhaled volume (FIF₇₅)

Forced expiratory flow at 50% of volume as a percentage of FIF50 (FEF50/FIF50)

The time taken between 25% and 75% of the forced expired volume (MET2575)

Forced Expiratory Time (FET)

Tidal Volume (TV)

Expiratory reserve volume (ERV)

Inspiratory reserve volume (IRV)

Inspiratory capacity (IC)

Expiratory Relaxed vital capacity (EVC)

Inspiratory vital capacity (IVC)

Breathing frequency rate (FR)

Inspiratory time (Ti)

Expiratory time (Te)

Ti as a % of total breath time (Ti/Ttot)

Tidal volume as a % of Ti (TV/Ti)

Breath Rate BR

Breathing Time B.T

Volume Tidal VT

Expiratory Time - average time of expiration per breath in seconds Te

Inspiratory Time - average time of inspiration per breath in seconds Ti

Total Tidal Breath Time in Seconds TTOT = Ti + Te

Ratio of Average Expiratory and Inspiratory Breaths Ti/Te

Average Time of Expiration per Breath as a ratio to the Total Tidal Breath Time Ti/TTOT

Tests per subject: 5 VC maneuver

8 FVC maneuvers

Predicted Values: Various – depends upon national preference

Transducer: CareFusion Bi-Directional Digital Volume.

Consumables/ Accessories

Cat. No.	Description
36-PSA1000	Adult Disposable Mouthpieces (500 per box)
36-SST1000	One-way Safety Mouthpieces (500 per box)
36-SST1250	One-way Safety Mouthpieces (250 per box)
36-PSA1200	Pediatric Disposable Mouthpieces (250 per box)
36-PSA1100	Pediatric Adaptor
V-892381	MicroGard® IIB with integrated mouthpiece (50 per box)
V-892384	MicroGard® IIC (50 per box)
36-SSC5000A	PeraSafe Sterilizing Powder 81g (to make up 5 liters of solution)
36-VOL2104	Nose Clips (pack of 5)
892121	Nose clip pads "foam material", 100 pcs.
36-SM2125	3 Liter Calibration Syringe
V-861449	Silicone adapter "oval" for MicroGard® IIB
861427	Calibration syringe adapter
36-SPC1000	Spirometry PC Software
36-CAB1047	USB Extension cable
36-ASS1141	SpiroUSB desktop stand
36-TDX1043	Turbine Transducer

For further information or to place an order for Disposables/
Supporting products please contact CAREFUSION, your local
distributor or view our website www.carefusion.com/micromedical

PLEASE NOTE: USE ONLY CAREFUSION ACCESSORIES

Please Note: Information in this manual is subject to change without notice and does not represent commitment on the part of CareFusion 232 UK Ltd. The software may be used or copied only in accordance with the terms of that agreement. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of CareFusion 232 UK Ltd.

Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551

Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.internationalsales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

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© CareFusion 2016
CareFusion Germany 234 GmbH