

DEMO, Ludo

PRE/POST ID


Masc. 53 Ans 95,0 kg 175,0 cm IMC: 31,02 Né(s) le: 29/03/1971

Valeurs prédites: C - GLI 2012; D - ECCS (Quanjer, 1993)


Capacité Vitale Forcée		PRE 28/01/2025 10:29:47					POST 28/01/2025 10:30:24					DIFFÉRENCE	
		53 ans, 175,0 cm, 95,0 kg, IMC : 31,0					53 ans, 175,0 cm, 95,0 kg, IMC : 31,0						
Résultats détaillés		Préd.	Min.	Mes.	%	Z-Sc.	Préd.	Min.	Mes.	%	Z-Sc.	Diff.	%
CVF	/ C	4,71	3,66	4,43	94	-0,43	4,71	3,66	4,38	93	-0,51	-0,05	-1
VEIMS	/ C	3,70	2,87	3,76	102	0,12	3,70	2,87	3,78	102	0,15	0,01	0
VEIMS/CVF	%	78,82	67,63	84,92	108	0,89	78,82	67,63	86,21	109	1,08	1,29	2
DEP	/s D	8,62	6,63	9,83	114	1,00	8,62	6,63	11,25	131	2,18	1,42	14
DEM75	/s D	7,55	4,74	8,38	111	0,49	7,55	4,74	9,40	124	1,08	1,01	12
DEM25-75	/s C	3,32	1,76	4,34	131	1,07	3,32	1,76	4,48	135	1,21	0,13	3
DEM50	/s D	4,64	2,47	4,70	101	0,05	4,64	2,47	4,62	100	-0,01	-0,08	-2
DEM25	/s D	1,85	0,57	1,42	77	-0,56	1,85	0,57	1,56	84	-0,38	0,14	10
TEF	s			6,28					5,92			-0,36	-6
AGE PULMO	ans			43,81					43,35			-0,46	-1

Capacité Vitale Lente		PRE 28/01/2025 10:13:07					POST					DIFFÉRENCE	
		53 ans, 175,0 cm, 95,0 kg, IMC : 31,0					53 ans, 175,0 cm, 95,0 kg, IMC : 31,0						
Résultats détaillés		Préd.	Min.	Mes.	%	Z-Sc.	Préd.	Min.	Mes.	%	Z-Sc.	Diff.	%
CV	/ C	5,11	4,09	4,42	86	-1,12							


CVF Courbe Débit/Volume




CVF Courbe Volume/Temps




CVL Courbe Volume/Temps



Miller Quadrant



Z-Sc. (PRE)



Z-Sc. (POST)

